



46241

Silver Proteinate extrapure, 8% Ag

Part B

CAS : 9008-42-8

Specifications

Appearance (Colour)	Brown
Appearance (Form)	Crystalline powder
Silver (Ag) content	~ 8.0%
Loss on drying	max. 7%

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	28432900 (GST 18%)
5 Gms	28432900 (GST 18%)

Available Packages

5 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com