



46448

Universal Restriction Enzyme Buffer (10X)

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Colourless
Appearance (Clarity)	Clear
Concentration	10X
Activity	Restriction digestion done with multiple restriction enzymes shows clear and complete digestion when incubated at 37°C for 1 hour.
Absence of nucleases	No degradation seen when incubated overnight with the buffer
Store at - 20°C	-20°C

Other Information

Description	An Unique buffer that can be used for most Restriction enzyme reaction set-ups. This buffer aides in optimal digestion within the specified time with great efficiency
Application	An Unique buffer that can be used for most Restriction enzyme reaction set-ups. This buffer aides in optimal digestion within the specified time with great efficiency

General Information

Storage	-20 °C (Blue/Dry Ice)
HSN Code	
5 x1 ml	35079099 (GST 18%)

Available Packages

5 x1 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com