



46821

Sodium Lauryl Sulphate (SLS, SDS) (Needle-shaped) extrapure, 90%

Part A

CAS : 151-21-3

Molecular Formula : C₁₂H₂₅SO₄Na

Molecular Weight : 288.38

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Needle shape
Solubility (Turbidity) 10% aq. solution	Clear to slight opalescent
Solubility (Colour) 10% aq. solution	Colourless to pale yellow
Assay (On dry basis)	min. 90%
pH (1% aq. solution)	8.0 -10.5
Alkalinity	max. 5ml N%
Sodium chloride + Sodium sulphate	max. 4%
Unsulphated alcohols	max. 2%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1325 , IMCO Class No.:4.1 , Packing Group:III
HSN Code	
5 Kg	34029019 (GST 18%)
500 Gms	34029019 (GST 18%)
Type of Packing	
5 Kg	Plastic Drum
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com