



46832

Amino Acids Kit

Part B

Specifications

Assay

All 24 Amino acids in the kit are >99% in purity

Homogeneity by TLC

All 24 Amino acids in the kit are chromatographically homogeneous by TLC

Other Information

Notes

22 amino acids - 1g each & 2 amino acids (L-Hydroxyproline, DL-DOPA) - 100mg each, in the kit, have assay min.99% and TLC test shows single spot.

Amino acids contained are,
DL-Alanine, DL-2-Aminobutyric acid,
L-Arginine monohydrochloride, DL-Aspartic acid,
L-Cysteine monohydrochloride, L-Cystine,
DL-DOPA, L-Glutamic acid,
Glycine, L-Histidine monohydrochloride,
L-Hydroxyproline, DL-Isoleucine, DL-Norleucine,
L-Leucine, L-Lysine monohydrochloride, DL-Methionine,
L-Ornithine monohydrochloride, DL-Phenylalanine,
L Proline, DL-Serine, DL-Threonine, L-Tryptophan,
L-Tyrosine, DL-Valine

General Information

Storage

25 to 40°C (Room Temperature)

Shelf Life

60 Months

IMDG Identification

Not Regulated for Transport (Non-Haz)

HSN Code

1 Kit

38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com