



47014

Arsenazo III Free Acid (High Purity) extrapure AR, 85%

Part B

CAS : 1668-00-4

Molecular Formula : C₂₂H₁₈As₂N₄O₁₄S₂

Molecular Weight : 776.37

Specifications

Appearance (Colour)	Brown to very dark brown and black
Appearance (Form)	Powder
Absorption maxima	569 - 576 nm
Absorbance (A) of 1% solution (pH 10.0) in a 1 cm cell	300 - 450
Assay (UV)	min. 85%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3465 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
1 Gms	29270090 (GST 18%)
5 Gms	29270090 (GST 18%)

Available Packages

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com