



47062

MES1: Spectrometric Oil Standard D12

Part B

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to Pale yellow
Appearance (Form)	Liquid
Silver (Ag)	10.00 ± 0.10 µg/g
Aluminium (Al)	10.00 ± 0.10 µg/g
Chromium (Cr)	10.00 ± 0.10 µg/g
Copper (Cu)	10.00 ± 0.10 µg/g
Iron (Fe)	10.00 ± 0.10 µg/g
Magnesium (Mg)	10.00 ± 0.10 µg/g
Sodium (Na)	10.00 ± 0.10 µg/g
Nickel (Ni)	10.00 ± 0.10 µg/g
Lead (Pb)	10.00 ± 0.10 µg/g
Silicon (Si)	10.00 ± 0.10 µg/g
Tin (Sn)	10.00 ± 0.10 µg/g
Titanium (Ti)	10.00 ± 0.10 µg/g

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
HSN Code	
100 Gms	38229090 (GST 12%)

Available Packages

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com