



47380

Polyoxyl 40 Hydrogenated Castor Oil extrapure

Part B

CAS : 61788-85-0

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Semisolid / liquid
Acid value	max 2.0
Saponification value	50 - 65mg KOH/g

Other Information

Application

This product is used to solubilize vitamins, water insoluble active substances and essential oils in water or mixtures of water and alcohol. Equivalent to Kolliphor RH40, Cremophor® RH 40.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	34024200 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com