



47448

## Antimony Potassium Tartrate Hemihydrate pure, 98%

Part A

CAS : 28300-74-5

Molecular Formula : C<sub>4</sub>H<sub>4</sub>O<sub>7</sub>KSb.1/2H<sub>2</sub>O

Molecular Weight : 333.93

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	min. 98%
Heavy Metals (Pb)	max. 0.002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1551 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
250 Gms	29181390 (GST 18%)
500 Gms	29181390 (GST 18%)
Type of Packing	
250 Gms	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

250 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com