



47471

Ethylene Glycol (MEG, Monoethylene Glycol)
pure, 98%

Part A

CAS : 107-21-1

Molecular Formula : C2H6O2

Molecular Weight : 62.07

Specifications	
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	1.112-1.114
Refractive Index (20°C)	1.431-1.432
Acidity	max. 0.1ml N%

General Information	
Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Ltr	29053100 (GST 18%)
2500 ml	29053100 (GST 18%)
500 ml	29053100 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages	
500 ml	
2500 ml	
25 Ltr	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com