



47484

Phenol Equilibrated with 0.1M Citrate Buffer pH 4.5 for molecular biology w/ Stabilizer

Part A

CAS : 108-95-2

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Yellow
Appearance (Form)	Liquid
pH	~ 4.5
DNase, RNase, Protease	Not detected
Stabilizer	Stabilized with 0.1% 8 -Hydroxy quinoline

Other Information

Notes	The product is suitable for RNA purification.
-------	---

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2821 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
100 ml	29071190 (GST 18%)
500 ml	29071190 (GST 18%)

Available Packages

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com