



47835

## Lanthanum (III) Trifluoromethanesulfonate extrapure, 98%

Part B

CAS : 52093-26-2

Molecular Formula : (CF<sub>3</sub>SO<sub>3</sub>)<sub>3</sub>La

Molecular Weight : 586.11

### Specifications

Appearance (Colour)	White to almost white
Appearance (Form)	Crystalline powder
Assay (Ion exchange)	min. 98%
Lanthanum content (La)	min. 23.2%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
1 Gms	28469090 (GST 18%)
5 Gms	28469090 (GST 18%)

### Available Packages

1 Gms  
5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com