



47994

## Paraffin Wax Pellets (Type 2 - 58-60) extrapure

Part A

CAS : 8002-74-2

### Specifications

Appearance (Colour)	White
Appearance (Form)	Soft pellets/flakes
Melting Point	58 - 60°C

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1000 Gms	27129090 (GST 18%)
500 Gms	27129090 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms  
1 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com