



48237

## DPPE (1,2-Bis(Diphenylphosphino)Ethane) extrapure, 98%

Part B

CAS : 1663-45-2

Molecular Formula : C<sub>26</sub>H<sub>24</sub>P<sub>2</sub>

Molecular Weight : 398.4

### Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Powder
Assay	min.98%
Melting Point	139 - 142°C

### Other Information

Application

Generally used in the following chemical reactions namely metal catalyzed allylic alkylation; decarboxylation of allylic esters, 1,3-diene synthesis, cycloaddition reactions, carbonylation reactions.

### General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	29031110 (GST 18%)
5 Gms	29031110 (GST 18%)

### Available Packages

5 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com