



48290

ProxiO 1kb DNA Ladder

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Blue
A260/A280	1.7~1.9
Concentration	100 - 120ng/ μ l
Stability	No obvious degradation was seen after incubation at 37°C for 1 week. Well-defined bands are formed during agarose gel electrophoresis.

Other Information

Description	The ProxiO 1kb ladder is suitable for sizing linear double-stranded DNA fragments from 500 bp-10 kb, and composed of ten chromatography-purified individual DNA fragments. All bands (except 2kb and 5kb) are supplied at approximately 40 ng/5mcl . The 2kb and 5kb bands have increased intensity relative to the other bands on an ethidium bromide-stained agarose gel.
Includes	50 mcg • ProxiO 1Kb Ladder - 50 mcg • 6X Orange Loading Dye - 1.0 ml 5 x50 mcg • ProxiO 1 Kb Ladder - 5 x 50 mcg • 6X Orange Loading Dye - 5 x 1.0 ml

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 x50 mcg	38229090 (GST 12%)
50 mcg	38229090 (GST 12%)

Available Packages

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com