



48656

## Cocktail 'T' scintillation grade

Part A

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pink fluorescent
Appearance (Form)	Solution
Absorbance of sample solution in a 1 cm cell v/s H <sub>2</sub> O @420nm	. max. 0.02

### Other Information

Notes	Content per litre : 5g PPO and 0.15g POPOP in toluene and Triton X-100
-------	---

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
1000 ml	38229090 (GST 12%)
2500 ml	38229090 (GST 12%)

### Available Packages

1000 ml
2500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com