



48914

Phenolphthalein Diphosphate Tetrasodium Salt extrapure AR, 95%

Part B

CAS : 68807-90-9

Molecular Formula : C₂₀H₁₂O₁₀P₂Na₄

Molecular Weight : 566.20

Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% aq. solution	Clear
Assay (NT) (On dry basis)	min. 95%
Free Phenolphthalein	max. 0.05%
Water	max. 16%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29322090 (GST 18%)
25 Gms	29322090 (GST 18%)
5 Gms	29322090 (GST 18%)

Available Packages

1 Gms

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com