



49092

## 2-Deoxycytidine extrapure, 98%

Part B

CAS : 951-77-9

Molecular Formula : C<sub>9</sub>H<sub>13</sub>N<sub>3</sub>O<sub>4</sub>

Molecular Weight : 227.22

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay (UV)	min. 98%
Spectral Ratio (A250nm/A260nm)	0.40 - 0.46
Spectral Ratio (A280nm/A260nm)	2.05 - 2.15
Spectral Ratio (A290nm/A260nm)	1.45 - 1.55

### General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29349990 (GST 18%)
100 Mg	29349990 (GST 18%)
500 Mg	29349990 (GST 18%)
Type of Packing	
1 Gms	Glass Bottle
100 Mg	Glass Bottle
500 Mg	Glass Bottle

### Available Packages

100 Mg  
500 Mg  
1 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com