



49244

Forchlorfenuron (KT-30, CPPU, N-(2-Chloro-4-pyridyl)-N-phenylurea) extrapure, 99%

Part B

CAS : 68157-60-8

Molecular Formula : C₁₂H₁₀ClN₃O

Molecular Weight : 247.68

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay (HPLC)	min. 99%
Melting Point	170 - 173°C
Loss on drying	max. 0.5%

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
1 Gms	98020000 (GST 18%)
5 Gms	98020000 (GST 18%)
Type of Packing	
1 Gms	Plastic Bottle
5 Gms	Glass Bottle

Available Packages

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com