



49456

## Vanillin extrapure AR, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 121-33-5

Molecular Formula : C<sub>8</sub>H<sub>8</sub>O<sub>3</sub>

Molecular Weight : 152.15

### Specifications

ExiPlus	meets compendial specs of BP
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in methanol	Clear
Solubility (Colour) 10% solution in methanol	Colourless
Assay (NT)	min. 99%
Melting Point	81- 83°C
Sulphated Ash	max. 0.01%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29124100 (GST 18%)
25 Gms	29124100 (GST 18%)
500 Gms	29124100 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle
500 Gms	Glass Bottle

### Available Packages

25 Gms

100 Gms

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com