



49473

Tetrasodium Pyrophosphate Decahydrate (TSPP Decahyd) extrapure AR, 99%

Part A

CAS : 13472-36-1

Molecular Formula : Na₄P₂O₇·10H₂O

Molecular Weight : 446.06

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay	min. 99%
pH (5% aq. solution)	9.5 - 10.5
Chloride (Cl)	max. 0.002%
Sulphate (SO ₄)	max. 0.005%
Iron (Fe)	max. 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2.5 Kg	28353900 (GST 18%)
500 Gms	28353900 (GST 18%)
Type of Packing	
2.5 Kg	Plastic Bottle
500 Gms	Glass Bottle

Available Packages

500 Gms

2.5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com