



49672

p-Xylenol Blue extrapure AR, 95%

Part A

CAS : 125-31-5

Molecular Formula : C₂₃H₂₂O₅S

Molecular Weight : 410.49

Specifications

Appearance (Colour)	Brown to Very Dark Brown
Appearance (Form)	Powder
Assay	min. 95%
Molar Extinction coefficient (?)	
@419-429nm	min. 8000
Visual transition interval	pH 2.8 Yellow

General Information

Storage	25 to 40°C (Room Temperature)
HSN Code	
1 Gms	32049000 (GST 18%)
5 Gms	32049000 (GST 18%)

Available Packages

1 Gms
5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com