



49982

Luminol extrapure, 98%

Part B

CAS : 521-31-3

Molecular Formula : C₈H₇N₃O₂

Molecular Weight : 177.16

Specifications

Appearance (Colour)	Yellow
Appearance (Form)	Coloured powder
Solubility (Turbidity) 1% solution in 1N NaOH	Clear
Assay	min. 98%
Loss on drying	max. 1%
Molar Extinction coefficient (?)	~7650
A300nm/A347nm	0.899 (±0.011)
IR Identification Test	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	29147990 (GST 18%)
5 Gms	29147990 (GST 18%)

Available Packages

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com