



50217

S-2-Chloropropionic Acid pure, 98%

Part B

CAS : 29617-66-1

Molecular Formula : CH₃CH(Cl)CO₂H

Molecular Weight : 108.5

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Yellow
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	1.257 - 1.262
Refractive Index (20°C)	1.432 - 1.434
Specific rotation [?] @20°C	-12° to -16°

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
100 ml	29155000 (GST 18%)
25 ml	29155000 (GST 18%)
5 ml	29155000 (GST 18%)
Type of Packing	
100 ml	Glass Bottle
5 ml	Glass Bottle

Available Packages

5 ml
25 ml
100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com