



50439

CdTe Core Type COOH Functionalized Quantum Dots, Fluorescence 510 nm (Hydrophilic)

Part C

Specifications		
Appearance (Colour)	Orange red	
Appearance (Form)	Compound	
Emission wavelength	510 nm	
General Information		
Storage	2 to 8°C (Refrigerate)	
Shelf Life	36 Months	
IMDG Identification	Not Regulated for Transport (Non-Haz)	
HSN Code		
10 Mg	28429020 (GST 18%)	
Type of Packing		
10 Mg	Glass Bottle	
Available Packages		
10 Mg		

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com