



50591

N,N-Dimethylformamide (DMF) for HPLC & UV Spectroscopy, 99.9%

Part A

CAS : 68-12-2

Molecular Formula : C₃H₇NO

Molecular Weight : 73.10

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99.9%
Density (g/ml) @ 20°C	0.947-0.949
Refractive Index (20°C)	1.429-1.431
Boiling Range	151-154°C
Non Volatile Matter	max. 0.0005%
Acidity (HCOOH)	max. 0.005%
Water (KF)	max. 0.05%
%T (Transmission) in a 1cm cell v/s H ₂ O	
@270nm	?30
@275nm	?65
@280nm	?75
@290nm	?85
@300nm	?90
@310nm	?95
@320nm	?97
@330nm	?98
@340-400nm	?99
(after purging with N ₂ for 15 min.)	

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2265 , IMCO Class No.:3 , Packing Group:III
HSN Code	
1000 ml	29211110 (GST 18%)

500 ml

29211110 (GST 18%)

Type of Packing

500 ml

Glass Bottle

Available Packages

500 ml

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
