



50615

## Tristearin pure, 90%(GC)

Part B

CAS : 555-43-1

Molecular Formula : C<sub>57</sub>H<sub>110</sub>O<sub>6</sub>

Molecular Weight : 891.51

### Specifications

Appearance (Colour)	White to off-white with yellow tinge
Appearance (Form)	Crystalline powder / flakes
Assay (GC)	min.90%
Melting Point	66 - 74°C

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29157020 (GST 18%)
25 Gms	29157020 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

### Available Packages

25 Gms  
100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com