



50716

## Iridium Trichloride Hydrate extrapure, 99%

Part A

CAS : 14996-61-3

Molecular Formula :  $\text{IrCl}_3 \cdot x\text{H}_2\text{O}$

Molecular Weight : 298.58 (anhydrous)

### Specifications

Appearance (Colour)	Grey to dark brownish black to dark Green
Appearance (Form)	Powder
Assay	min. 99%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3260 , IMCO Class No.:8 , Packing Group:II
HSN Code	
1 Gms	28439019 (GST 18%)
250 Mg	28439019 (GST 18%)
5 Gms	28439019 (GST 18%)

### Available Packages

250 Mg

1 Gms

5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com