



51113

Inosine extrapure, 99%

Part B

CAS : 58-63-9

Molecular Formula : C₁₀H₁₂N₄O₅

Molecular Weight : 268.23

Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Assay (UV)	min. 99%
Heavy Metals (Pb)	max. 0.001%
Spectral Ratio (A250nm/A260nm)	1.73 ± 0.10
Spectral Ratio (A280nm/A260nm)	0.24 ± 0.05

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29349990 (GST 18%)
100 Gms	29349990 (GST 18%)
25 Gms	29349990 (GST 18%)
500 Gms	29349990 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

10 Gms
25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com