



51126

Chloroform Dried, 99.5%, water 0.005%

Part A

CAS : 67-66-3

Molecular Formula : CHCl₃

Molecular Weight : 119.38

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	1.475 - 1.479
Refractive Index (20°C)	1.444 - 1.445
Boiling Range	60 - 62°C
Non Volatile Matter	max. 0.0005%
Acidity (HCl)	max. 0.0005%
Aldehydes & Ketones (CH ₃ CHO)	max. 0.005%
Substance Discoloured by H ₂ SO ₄	Passes test
Organic impurities	Passes test
Dithizone test	Passes test
Metallic impurities	No reaction
Ionised Chloride (Cl)	max. 0.00002%
Free Chlorine (Cl)	max. 0.00002%
Iron (Fe)	max. 0.00001%
Heavy Metals (Pb)	max. 0.00001%
Water (KF)	max. 0.005%
Stabilizer (Amylene)	Present

Other Information

Notes

Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1888 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
1000 ml	29031300 (GST 18%)
2500 ml	29031300 (GST 18%)
500 ml	29031300 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle

Available Packages

500 ml
1000 ml
2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
