



51171

4-Carboxyphenylboronic Acid extrapure, 97%

Part B
CAS : 14047-29-1
Molecular Formula : C7H7BO4
Molecular Weight : 165.9

Specifications	
Appearance (Colour)	White to off-white
Appearance (Form)	Powder
Assay	min. 97%
Melting Point	220 - 225°C

General Information	
Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
HSN Code	
1 Gms	29319090 (GST 18%)
5 Gms	29319090 (GST 18%)
Type of Packing	
1 Gms	Plastic Bottle
5 Gms	Glass Bottle

Available Packages	
1 Gms	
5 Gms	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.
608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com