



51451

Bismuth (III) Chloride pure, 99%

Part A

CAS : 7787-60-2

Molecular Formula : BiCl3

Molecular Weight : 315.34

Specifications

Appearance (Colour)	White to biege to gray
Appearance (Form)	Powder
Assay (Complexometric)	min. 99%
Sulphate (SO4)	max. 0.02%
Iron (Fe)	max. 0.02%
Nitrate (NO3)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3260 , IMCO Class No.:8 , Packing Group:III
HSN Code	
100 Gms	28273990 (GST 18%)

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com