



51535

Silver Nitrate 0.1N Solution

Part A

CAS : 7761-88-8

Molecular Formula : AgNO₃

Molecular Weight : 169.87

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Normality (N)	0.1000 ± 0.005

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	36 Months
IMDG Identification	UN No.:1493 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
1000 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)
Type of Packing	
1000 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml

1000 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com