



51775

Ammonium Iodide extrapure, 99%

Part A

CAS : 12027-06-4

Molecular Formula : NH₄I

Molecular Weight : 144.94

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Crystalline powder
Assay (T)	min. 99.0%
Residue on ignition	max. 0.1%
Sulphate (SO ₄)	max. 0.01%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28276090 (GST 18%)
250 Gms	28276090 (GST 18%)
50 Gms	28276090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
250 Gms	Glass Bottle
50 Gms	Glass Bottle

Available Packages

50 Gms
100 Gms
250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com