



51849

Silica Gel G for TLC, 325mesh, w/ CaSO4

Part A

CAS : 112926-00-8

Molecular Formula : SiO₂

Specifications

Appearance (Colour)	White
Appearance (Form)	Free flowing powder
pH (10% aq. solution)	~ 7.0
Chloride (Cl)	max. 0.02%
Iron (Fe)	max. 0.02%
Heavy Metals (Pb)	max. 0.02%
Suitability for TLC	Passes test
Adhesive Property	Passes test
Organic impurities	Passes test
Binder (CaSO ₄ ·1/2H ₂ O)	~ 13%
Pore size	60°A

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Kg	28112200 (GST 18%)
5 Kg	28112200 (GST 18%)
500 Gms	28112200 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms
5 Kg
25 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com