



51952

## Pyridine Dried, 99.5%, water 0.01%

Part A

CAS : 110-86-1

Molecular Formula : C<sub>5</sub>H<sub>5</sub>N

Molecular Weight : 79.10

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	0.980-0.983
Refractive Index (20°C)	1.508-1.510
Boiling Range	114-116°C
Non Volatile Matter	max. 0.002%
Water (KF)	max. 0.01%

### Other Information

Notes

Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1282 , IMCO Class No.:3 , Packing Group:II
HSN Code	
500 ml	29333100 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

### Available Packages

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com