



52208

FMOC-O-Benzyl-L-Serine extrapure, 98%

Part B

CAS : 83792-48-7

Molecular Formula : C₂₅H₂₃NO₅

Molecular Weight : 417.45

Specifications

Appearance (Colour)	White to Light Yellow
Appearance (Form)	Powder to crystals
Solubility (Turbidity) 1% solution in methanol	Clear
Assay (NaOH titration)	min. 98%
Specific rotation [α] _D ²⁰	+23° to +27 (C=1, ethyl acetate)
Melting Point	144-148°C
Water (KF)	max.1.0%

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)

Available Packages

10 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com