



52519

Boron Trifluoride Diethyl Etherate pure, 45-50% BF3

Part A

CAS : 109-63-7

Molecular Formula : BF₃.C₄H₁₀O

Molecular Weight : 141.93

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow
Appearance (Form)	Liquid
Assay (T)	min. 45 - 50% BF ₃
Boiling Range	126-129°C

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2604 , IMCO Class No.:8(3) , Packing Group:I
HSN Code	
2500 ml	29420090 (GST 18%)
500 ml	29420090 (GST 18%)

Available Packages

250 ml

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com