



52693

Sodium Metabisulphite extrapure AR, ACS, ExiPlus, Multi-Compendial, 98%

Part A

CAS : 7681-57-4

Molecular Formula : Na₂S₂O₅

Molecular Weight : 190.10

Specifications

ExiPlus	meets compendial specs of USP, BP
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay	min. 98%
pH (5% aq. solution)	3.8 - 5.3
Chloride (Cl)	max. 0.005%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.001%
Arsenic (As)	max. 0.0001%
Thiosulphate (S ₂ O ₃)	max. 0.02%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	28321090 (GST 18%)
500 Gms	28321090 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com