



53032

## Tungsten Disulfide Micronpowder, 0.6micron, 99.9%

Part C

CAS : 12138-09-9

Molecular Formula : WS<sub>2</sub>

Molecular Weight : 247.98

### Specifications

Appearance (Colour)	Grey to black
Appearance (Form)	Powder
Assay	min. 99.9%
APS	0.6 micron

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28309020 (GST 18%)
25 Gms	28309020 (GST 18%)
500 Gms	28309020 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

25 Gms  
100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com