



53061

FMOC-L-Amino Acids Kit

Part B

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	All FMOC-L-Amino acids in the kit are min 99%

Other Information

Notes

All 20 FMOC-L-Amino acids, 5g each, in the kit, have assay min.99% and TLC test shows single spot.

Amino acids contained are,
FMOC-L-Alanine, FMOC-L-Arginine,
FMOC-L-Asparagine, FMOC-L-Citrulline,
FMOC-L-Glutamine, FMOC-L-Glutamic
Acid,
FMOC-L-Histidine, FMOC-L-
Isoleucine,
FMOC-L-Leucine, FMOC-L-Lysine,
FMOC-L-Methionine, FMOC-L-
Norleucine,
FMOC-L-Norvaline, FMOC-L-
Phenylalanine,
FMOC-L-Proline, FMOC-L-Serine,
FMOC-L- Threonine, FMOC-L-Trypt

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
HSN Code	
1 Kit	38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com