



53110

D-Lactate Dehydrogenase (D-LDH) ex. Microorganism, 400U/mg solids

Part B

CAS : 9028-36-8

Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Activity	min. 400 units/mg powder
Unit Definition	One unit will reduce 1.0 micromole of pyruvate to D-lactate per min at pH 7.0 at 25 °C.
Enzyme solution for activity determination	Prepared in 0.1 M potassium phosphate buffer, pH 7.0.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2500 Units	35079099 (GST 18%)
5000 Units	35079099 (GST 18%)

Available Packages

2500 Units

5000 Units

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com