



53180

1,1,3,3-Tetramethylguanidine extrapure AR, 99%

Part B

CAS : 80-70-6

Molecular Formula : C5H13N3

Molecular Weight : 115.18

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.916-0.920
Refractive Index (20°C)	1.467-1.469

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
HSN Code	
100 ml	29252990 (GST 18%)
500 ml	29252990 (GST 18%)

Available Packages

100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com