



53422

## Diethyl Ketone (3-Pentanone) extrapure, 99%

Part A

CAS : 96-22-0

Molecular Formula : C<sub>5</sub>H<sub>10</sub>O

Molecular Weight : 86.13

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.813-0.815
Refractive Index (20°C)	1.391-1.392
Boiling Range	101-102°C
Acidity (CH <sub>3</sub> COOH)	max. 0.05%
Water (KF)	max. 0.2%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1156 , IMCO Class No.:3 , Packing Group:II
HSN Code	
500 ml	29141990 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

### Available Packages

500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com