



53457

D-Serine extrapure, 99%

Part A

CAS : 312-84-5

Molecular Formula : C₃H₇NO₃

Molecular Weight : 105.1

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% solution 2N HCl	Clear
Solubility (Colour) 10% solution 2N HCl	Colourless
Assay	min. 99%
Specific rotation [?]	-13° to -14° (C=5 , 5N HCl)
Loss on drying	max. 0.3%
Sulphated Ash	max. 0.1%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%
Homogeneity by TLC	Passes test
Not for medicinal use	

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29224990 (GST 18%)
25 Gms	29224990 (GST 18%)
5 Gms	29224990 (GST 18%)
500 Gms	29224990 (GST 18%)

Available Packages

25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com