



53534

FMOC-D-Serine extrapure, 99%

Part B

CAS : 116861-26-8

Molecular Formula : C₁₈H₁₇NO₅

Molecular Weight : 327.34

Specifications

Appearance (Colour)	White to slight yellow to beige
Appearance (Form)	Powder
Solubility (Turbidity) 10% solution in DMF	Clear
Solubility (Colour) 10% solution in DMF	Colourless to pale yellow
Assay (NT) (On dry basis)	min 99.0%
Water (KF)	max 7.0%

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2.5 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
500 Mg	29379020 (GST 18%)

Available Packages

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com