



53625

ProxiB Lambda Double Digest Marker : 0.5-23kb

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Blue
Absence of nucleases	No degradation seen after incubation at 37°C for 24 hours
Agarose gel electrophoresis	Well defined bands with clear intensity
Concentration	~100ng/μl

Other Information

Description	Lambda DNA-digested with HindIII and EcoRI is a ready-to load preparation and supplied as 1 mcg DNA / 10 mcl . HindIII, EcoRI double digestion of lambda DNA generates 12 fragments (in kb) 21.2, 5.1, 4.9, 4.2, 3.5, 2.0, 1.9, 1.5, 1.3, 0.9, 0.8, and 0.5.
-------------	--

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 mcg	38229090 (GST 12%)
50 mcg	38229090 (GST 12%)

Available Packages

50 mcg

100 mcg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com