



53713

## Sodium Biselenite pure, 98%

Part B

CAS : 7782-82-3

Molecular Formula : NaHSeO3

Molecular Weight : 150.96

### Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Assay (Iodometric)	min. 98%
pH (2% aq. solution)	4.0 - 7.5

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2630 , IMCO Class No.:6.1 , Packing Group:I
HSN Code	
100 Gms	28419000 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

### Available Packages

100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com