



53755

Digitonin (Digitin) extrapure, 50%

Part B

CAS : 11024-24-1

Molecular Formula : C₅₆H₉₂O₂₉

Molecular Weight : 1229.30

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 2.5% solution in Methanol	Clear
Solubility (Colour) 2.5% solution in Methanol	Colourless
Melting Point	230 - 240°C (Decomposition)
Specific rotation [?]	- 47° to -49° (C=1 ,75% Acetic acid)

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:3462 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
1 Gms	29420090 (GST 18%)
100 Mg	29420090 (GST 18%)

Available Packages

100 Mg

1 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com