



54236

Progoitrin Potassium Salt Reference Standard
Grade, 75%

Part B

CAS : 21087-77-4

Molecular Formula : C11H18KNO10S2

Molecular Weight : 427.48

Specifications

Appearance (Colour)	Off - white
Appearance (Form)	Hygroscopic powder
Assay	min. 75%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
10 Mg	29420090 (GST 18%)

Available Packages

10 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com