



54422

Poly-L-Lysine Hydrobromide (PLL HBr), 70-150kDa

Part B

CAS : 25988-63-0

Molecular Weight : 70,000-150,000

Molecular Formula : L-Lys-(L-Lys)_n-L-Lys.xHBr

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Solid
Solubility (Colour) 5% aq. solution	Colourless to faint yellow
Solubility (Turbidity) 5% aq. solution	Clear to slightly hazy
Identification (IR Spectrum)	Conforms to structure
Water (KF)	max.10%
Degree of Polymerization	330 - 720 (By 5% aq. viscosity)
Molecular weight by Viscosity	70 - 150KDa

Available Packages

25 Mg

100 Mg

250 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com